

Paraffin Dispenser XH

- The full automatic procedure control & new heating element gives feature of quick heating & energy saving
- The temperature survey integration block from American DALLAS Corporation has been used to make it precise in temperature survey and credible in performance
- It possesses the function of memorization to keep the set temperature automatically after startup
- The LED shows heating state



Technical Specifications

- Range of temperature: ambient~100°C
- Precision: ±1%
- Voltage: AC220V 50HZ / AC110V 60HZ
- Power: 1200VA
- Capacity: 10000ml
- Overall dimension: 350x450x510mm
- Net Weight: 17kg