



[Vented Storage Cabinets]

Baygen Vented Tissue Storage Cabinet is designed to connect ventilation system. The perforated shelves allows to removes formalin fumes. Thanks to this cabinet personnel are protected as harmful formalin fumes continually exhausted from the cabinets.



- Organized and Odorless Storage
- Short or Long Term Storage
- Connects to the Existing Ventilation System
- With solid or sliding doors



Technical Features

- All stainless steel
- AISI 304 stainless steel shelves and waste liquid collector
- Durable shelf system up to 100 kg
- Perforated shelves suitable for effective ventilation with 1.5cm diameter holes on 7.5 centers
- 16x5cm ventilation cuff is centered on the top of the cabinet
 - eliminate odors
 - short and long term specimen storage
 - be integrated to central ventilation system
- Solid or sliding doors

Models: 01-0005-01 Vented Storage Cabinet without Exhaust Fan
 01-0005-02 Vented Storage Cabinet with Exhaust Fan
 01-0005-03 Vented Storage Cabinet, Ductless

